

Code No.08-0910

# SPECIFICATION

## Hexane, Super Dehydrated for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.001%
Residue after evaporation	max.0.001%
Acidity (as CH <sub>3</sub> COOH)	max.0.002%
Fat	to pass test
Thiophene group	to pass test
Substances darkened by sulfuric acid	to pass test
Sulfide and reducing substances	to pass test
Assay (capillary GC)	min.96.0%