

Code No.08-1045

## SPECIFICATION

### Hyaluronic Acid Sodium Salt, from Rooster Comb for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellow, crystalline powder ~
Solubility in water	powder or granules to pass test
Loss on drying at 60°C over phosphorus(V) oxide in vacuum	max.15.0%
Residue after ignition (as sulfate) (after drying)	15.0~20.0%
pH of a 1g/L solution at 25°C	5.0~8.5
Nitrogen compounds (as N)	2.8~4.0%