

Code No.09-0311

SPECIFICATION

2,2'-Iminodiethanol

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [HN(CH ₂ CH ₂ OH) ₂] (GC) Solubility in water Refractive index 20°C Water Residue after ignition 2-Aminoethanol (H ₂ NCH ₂ CH ₂ OH) (GC) 2,2',2''-Nitrilotriethanol [N(CH ₂ CH ₂ OH) ₃] (GC)	Colorless ~ nearly colorless, clear liquid min.99.0%(mass/mass) to pass test 1.474~1.479 max.0.3%(mass/mass) max.0.01%(mass/mass) max.0.3%(mass/mass) max.0.3%(mass/mass)

JIS K 8453:2013(*:Additional test performed by Wako but not required by JIS)