

Code No.09-0311

SPECIFICATION

2,2'-Iminodiethanol

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ nearly colorless, clear liquid
Assay [HN(CH ₂ CH ₂ OH) ₂] (GC)	min.99.0%(mass/mass)
Solubility in water	to pass test
Refractive index 20°C	1.474~1.479
Water	max.0.3%(mass/mass)
Residue after ignition	max.0.01%(mass/mass)
2-Aminoethanol (H ₂ NCH ₂ CH ₂ OH) (GC)	max.0.3%(mass/mass)
2,2',2''-Nitrilotriethanol [N(CH ₂ CH ₂ OH) ₃] (GC)	max.0.3%(mass/mass)

JIS K 8453:2013(*:Additional test performed by Wako but not required by JIS)