

Code No.09-0524

SPECIFICATION

Interferon-lambda 1, Human, recombinant  
for Cellbiology

| REQUIREMENT                     | SPECIFICATION      |
|---------------------------------|--------------------|
| Appearance                      | to pass test       |
| Bacterial endotoxins            | max.0.1ng/ $\mu$ g |
| Electrophoresis test (SDS-PAGE) | to pass test       |
| Mass spectrometry (MALDI-TOF)   | to pass test       |

Source:E.coli expressed human interferon-lambda 1