

Code No.09-0698

## SPECIFICATION

### Interferon- $\gamma$ , Mouse, recombinant, Animal-derived-free for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ $\mu$ g
Electrophoresis test (SDS-PAGE)	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source: E. coli expressed mouse interferon-  $\gamma$

This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.