

SPECIFICATION

L-Leucine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₆ H ₁₃ NO ₂) Solubility in dil.hydrochloric acid pH of a 10g/L solution at 25°C Specific rotation [α] _{D20} Loss on drying at 105°C Residue after ignition (as sulfate) Chloride (Cl) Sulfate (SO ₄) Lead (Pb) Arsenic (As) Iron (Fe) Ammonium (NH ₄) Foreign amino acids	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 5.5~7.0 +14.7~+16.0° max.0.2%(mass/mass) max.0.05%(mass/mass) max.0.01%(mass/mass) max.0.02%(mass/mass) max.0.001%(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass) max.0.02%(mass/mass) to pass test

JIS K 9054 : 2020(*:Additional test performed by Wako but not required by JIS)