

## SPECIFICATION

## Liquid Paraffin

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Density (20°C)	0.86~0.89g/ml
Acidity (as H <sub>2</sub> SO <sub>4</sub> )	max.0.002%(mass/mass)
Alkalinity (as NaOH)	max.0.001%(mass/mass)
Chloride (Cl)	max.0.001%(mass/mass)
Sulfate (SO <sub>4</sub> )	max.0.001%(mass/mass)
Sulfide	to pass test
Solid paraffin	to pass test
Substances darkened by sulfuric acid	to pass test
Substances reducing permanganate (as O)	to pass test

JIS K 9003:2014

(\*:Additional test performed by Wako but not required by JIS)