

Code No.12-0532

SPECIFICATION

Ligroin

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Boiling range (80~110°C)	min.90%(volume/volume)
Residue after evaporation	max.0.001%(mass/mass)
Acidity (as H ₂ SO ₄)	max.5ppm(mass/mass)
Sulfide and reducing substances	to pass test
Substances darkened by sulfuric acid	to pass test
Copper corrosion test	to pass test

JIS K 8937:2020(*:Additional test performed by Wako but not required by JIS)