

# SPECIFICATION

## Loganin

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
Assay (qNMR) (JP method)	under investigation
* Solubility in methanol	to pass test
* Loss on drying over silica gel	max.1%

This product is suited for Loganin for assay (for assay 2) and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (\*:Additional test performed by Wako)