

## SPECIFICATION

## Magnesium Acetate Tetrahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder or mass
Solubility in dil.acetic acid	to pass test
Chloride (Cl)	max.0.002%
Phosphate (PO4)	max.5ppm
Sulfate (SO4)	max.0.01%
Nitrogen compounds (as N)	max.0.003%
Sodium (Na)	max.0.01%
Potassium (K)	max.0.005%
Copper (Cu)	max.5ppm
Calcium (Ca)	max.0.01%
Barium (Ba)	max.0.005%
Lead (Pb)	max.0.001%
Iron (Fe)	max.5ppm
Assay	min.99.0%