

SPECIFICATION

Manganese(II) Chloride Tetrahydrate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--|
| * Appearance Assay (MnCl ₂ ·4H ₂ O) Solubility in water pH of a 50g/L solution at 25°C Sulfate (SO ₄) Sodium (Na) Potassium (K) Copper (Cu) Magnesium (Mg) Calcium (Ca) Zinc (Zn) Lead (Pb) Iron (Fe) Nickel (Ni) | Pale red, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 4.0~6.5 max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass) max.5ppm(mass/mass) max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass) |

JIS K 8160:2017

(*:Additional test performed by Wako but not required by JIS)