

SPECIFICATION

Dichloromethane for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 300

| REQUIREMENT | SPECIFICATION |
|---|--|
| * Appearance Assay (CH ₂ Cl ₂) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis * Suitability for pesticide residue analysis (GC-FPD) * Methanol (stabilizer) | Colorless clear liquid min.99.5%(mass/mass) 1.315~1.330g/mL max.0.01%(mass/mass) max.0.0002%(mass/mass) to pass test to pass test 0.2~0.5%(mass/mass) |

JIS K 8117:2015(*:Additional test performed by Wako but not required by JIS)