

SPECIFICATION

3-Methyl-1-butanol

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [(CH ₃) ₂ CHCH ₂ CH ₂ OH] (GC) Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as CH ₃ COOH) Aldehydes Basic substances (as N) Substances darkened by sulfuric acid	Colorless clear liquid min.98.0%(mass/mass) 0.808~0.813g/ml 1.402~1.407 max.0.3%(mass/mass) max.0.002%(mass/mass) max.0.01%(mass/mass) to pass test max.0.001%(mass/mass) to pass test

JIS K 8051:2010(*:Additional test performed by Wako but not required by JIS)