

# SPECIFICATION

## Multielement Standard Solution W-V for ICP Analysis

REQUIREMENT	SPECIFICATION
Appearance	Pale blue clear liquid
Acidity (HNO <sub>3</sub> )	0.90~1.10mol/l
Aluminium (Al)	95~105mg/l
Boron (B)	95~105mg/l
Cadmium (Cd)	95~105mg/l
Chromium (Cr)	95~105mg/l
Copper (Cu)	95~105mg/l
Iron (Fe)	95~105mg/l
Manganese (Mn)	95~105mg/l
Molybdenum (Mo)	95~105mg/l
Sodium (Na)	95~105mg/l
Nickel (Ni)	95~105mg/l
Lead (Pb)	95~105mg/l
Zinc (Zn)	95~105mg/l