

SPECIFICATION

Multielement Standard Solution W-VI

for ICP Analysis

REQUIREMENT	SPECIFICATION
Appearance	Pale blue clear liquid
Acidity (HNO ₃)	0.09~0.11mol/l
Aluminium (Al)	95.0~105.0mg/l
Boron (B)	95.0~105.0mg/l
Calcium (Ca)	95.0~105.0mg/l
Cadmium (Cd)	95.0~105.0mg/l
Chromium (Cr)	95.0~105.0mg/l
Copper (Cu)	95.0~105.0mg/l
Iron (Fe)	95.0~105.0mg/l
Magnesium (Mg)	95.0~105.0mg/l
Manganese (Mn)	95.0~105.0mg/l
Sodium (Na)	95.0~105.0mg/l
Lead (Pb)	95.0~105.0mg/l
Zinc (Zn)	95.0~105.0mg/l