

SPECIFICATION

Multication Standard Solution I for Ion Chromatography

| REQUIREMENT | SPECIFICATION |
|-----------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Acidity (as HNO ₃) | 0.018~0.022mol/l |
| Sodium ion (Na ⁺) | 95.0~105.0mg/l |
| Magnesium ion (Mg ²⁺) | 95.0~105.0mg/l |
| Calcium ion (Ca ²⁺) | 95.0~105.0mg/l |