

SPECIFICATION

Multianion Standard Solution III for Ion Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Bromide ion (Br ⁻)	47.5~52.5mg/L
Chloride ion (Cl ⁻)	9.5~10.5mg/L
Nitrite ion (NO ₂ ⁻)	47.5~52.5mg/L
Nitrate ion (NO ₃ ⁻)	47.5~52.5mg/L
Sulfate ion (SO ₄ ²⁻)	95~105mg/L