

Code No.13-1534

SPECIFICATION

2-Methyl-1-propanol

for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	0.800~0.804g/mL
Refractive index 20°C	1.395~1.398
Water	max.0.2%
Acidity [as (CH ₃) ₂ CHCOOH]	max.0.002%
Fluorescence test	to pass test
Substances darkened by sulfuric acid	to pass test
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Assay (GC)	min.99.0%