

Code No.13-1620

SPECIFICATION

Multielement Standard Solution W-X

for ICP-MS Analysis

REQUIREMENT	SPECIFICATION
Appearance	Slightly blue, clear liquid
Acidity (HNO ₃)	0.18~0.22mol/l
Aluminium (Al)	95.0~105.0mg/l
Arsenic (As)	95.0~105.0mg/l
Cadmium (Cd)	95.0~105.0mg/l
Chromium (Cr)	95.0~105.0mg/l
Copper (Cu)	95.0~105.0mg/l
Manganese (Mn)	95.0~105.0mg/l
Lead (Pb)	95.0~105.0mg/l
Selenium (Se)	95.0~105.0mg/l
Zinc (Zn)	95.0~105.0mg/l