

## SPECIFICATION

## Magnolol

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
* Assay (qNMR) (JP method)	min.98.0%
** Solubility in methanol	to pass test
** Loss on drying over silica gel	max.1.0%

This product is suited for Magnolol for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia.  
(\*:The specification is originally set by Wako)(\*\*):Additional test performed by Wako)