

Code No.13-1939

SPECIFICATION

100mmol/L N1-Methylpseudouridine 5'-Triphosphate Sodium Solution  
for Nucleic Acid Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Liquid
pH	7.0~9.0
Calcium (Ca)	under investigation
Magnesium (Mg)	under investigation
Bacterial endotoxins	less than 1EU/mL
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Concentration	90~110mmol/L
Assay (HPLC)	min.98.5%