

Code No.14-0234

## SPECIFICATION

## Disodium p-Nitrophenylphosphate Hexahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	8.0~9.0
Chloride (Cl)	max.0.005%
Phosphate (PO <sub>4</sub> )	max.0.04%
Heavy metals (as Pb)	max.0.003%
p-Nitrophenol (O <sub>2</sub> NC <sub>6</sub> H <sub>4</sub> OH)	max.0.01%
Assay	min.97.0%