

## SPECIFICATION

6 Nitrosamines Mixture Standard Solution (each 2  $\mu\text{g}/\text{mL}$  Methanol Solution)  
for Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Practical test	to pass test
N-Nitrosodiethanolamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosopyrrolidine	1.97~2.03 $\mu\text{g}/\text{mL}$
4-(Methylnitrosoamino)-1-(3-pyridinyl)-1-butanone	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosopiperidine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Methyl-N-nitrosophenethylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodiphenylamine	1.97~2.03 $\mu\text{g}/\text{mL}$