

Code No.15-0220

SPECIFICATION

OxoneR Monopersulfate Compound

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystalline powder ~
Solubility in water	powder or granules to pass test
pH of a 50g/l solution at 25°C	1.4~1.9
Consumption of 0.1mol/l sodium hydroxide solution	11.0~21.0ml/g
Magnesium carbonate (MgCO ₃)	0.5~2.5%
Available oxygen (O)	4.5~5.1%