

Code No.15-0363

## SPECIFICATION

Oxidizing Solution [Iodine Solution (abt.  
0.02mol/L)][Tetrahydrofuran:Pyridine:Water(78:20:2)]  
for Nucleic Acid Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Brown liquid
Water	1.8~2.2vol%
Iodine (as I <sub>2</sub> )	0.015~0.024mol/L
Pyridine	19.0~21.0vol%
Tetrahydrofuran	76.0~80.0vol%