

Code No.16-0055

# SPECIFICATION

## Pectin, from Citrus

### Practical Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale brown, powder
Loss on drying at 105°C	max.15.0%
Residue after ignition (as sulfate) (after drying)	max.5.0%
Heavy metals (as Pb)	max.0.01%
Arsenic (as As)	max.2ppm
Methoxy group (CH <sub>3</sub> O-) (calculated on the dried basis)	3.5~8.5%
Galacturonic acid(C <sub>6</sub> H <sub>10</sub> O <sub>7</sub> ) (calculated on the dried basis)	50.0~80.0%