

## SPECIFICATION

## Potassium Acetate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH <sub>3</sub> COOK) Solubility in water pH of a 50g/l solution at 25°C Loss on drying at 150°C Chloride (Cl) Phosphate (PO <sub>4</sub> ) Sulfate (SO <sub>4</sub> ) Sodium (Na) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe) Substances reducing permanganate (as O)	White, crystals ~ crystalline powder min.97.0%(mass/mass) to pass test 7.5~9.0 max.3.0%(mass/mass) max.0.003%(mass/mass) max.0.001%(mass/mass) max.0.003%(mass/mass) max.0.5%(mass/mass) max.5ppm(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass)

JIS K 8363:2011

(\*:Additional test performed by Wako but not required by JIS)