

SPECIFICATION

Potassium Hexacyanoferrate(II) Trihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [K ₄ [Fe(CN) ₆]·3H ₂ O] Solubility in water Chloride (Cl) Sulfate (SO ₄) Sodium (Na)	Pale yellowish brown, crystals ~ crystalline powder min.99.5%(mass/mass) to pass test max.0.003%(mass/mass) max.0.005%(mass/mass) max.0.01%(mass/mass)

JIS K 8802:2022(*:Additional test performed by Wako but not required by JIS)