

Code No.16-0908

## SPECIFICATION

## Polyethylene Glycol 1,000

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, mass
Viscosity (500g/l, 25°C)	15~35mPa·s
Solubility in water	to pass test
Freezing point	30~40°C
Water	max.2.0%
Residue after ignition (as sulfate)	max.0.2%
pH of a 50g/l solution at 25°C	4.0~7.0
Acidity (as CH <sub>3</sub> COOH)	max.0.1%
Average molecular weight	900~1,100