

SPECIFICATION

Potassium Dihydrogenphosphate for Oxidant Analysis
for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance Assay (KH ₂ PO ₄) Solubility in water pH of a 50g/L solution at 25°C Nitrogen compounds (as N) Lead (Pb) Suitability for oxidant determination	White crystals min.99.0%(mass/mass) to pass test 4.2~4.6 max.0.001%(mass/mass) max.0.001%(mass/mass) to pass test

JIS K 8118:2021(*:Additional test performed by Wako but not required by JIS)