

# SPECIFICATION

## Potassium Gluconate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystalline powder or shot
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	under investigation
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.01%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Reducing sugars (as C <sub>6</sub> H <sub>12</sub> O <sub>6</sub> )	max.1.0%
Assay	min.97.0%