

## SPECIFICATION

Palladium(II) Nitrate Solution  
for Atomic Absorption Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Brown clear liquid
Acidity (HNO <sub>3</sub> )	under investigation
Silver (Ag)	max.50ppb
Cadmium (Cd)	max.1ppb
Thallium (Tl)	max.10ppb
Lead (Pb)	max.10ppb
Arsenic (As)	max.10ppb
Selenium (Se)	max.100ppb
Concentration (Pd) (20°C)	4.95~5.05mg/mL