

Code No.16-1925

## SPECIFICATION

### 2% Potassium Hydroxide-impregnated Silica Gel for Dioxins Analysis

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (75~150 $\mu$ m)	min.75%
Particle size (pass 45 $\mu$ m)	max.5%
Precipitation volume	1.6~2.3ml/g
Loss on drying at 105°C	max.2.0%
Practical test	to pass test
Potassium hydroxide (KOH) (calculated on the dried basis)	1.6~2.4%