

Code No.16-2191

SPECIFICATION

Prednisolone Standard

for High Performance Liquid Chromatography

| REQUIREMENT | SPECIFICATION | |
|---|--|------------|
| Appearance | White ~ nearly white, crystalline powder ~ powder to pass test | |
| Solubility in ethanol | | |
| Specific rotation [α]D20(c=1,C2H5OH) (after drying) | | +113~+119° |
| Loss on drying at 105°C | | max.1.0% |
| Assay (HPLC) | | min.98.0% |