

Code No.16-2734

SPECIFICATION

5% Phenyl Isothiocyanate, n-Heptane Solution

for SHIMADZU Protein Peptide Sequencer

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Diphenylthiourea	max.0.0005%
Phenylthiourea	max.0.0005%
Peroxides	to pass test
Suitability for amino acids sequence analysis	to pass test
Purity(HPLC)(Phenyl isothiocyanate)	min.99.9%
Assay	5.0~5.5%