

Code No.18-0268

# SPECIFICATION

## Resibufogenin

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystalline powder
Absorbance E1%1cm (300nm) (JP method)	131~145
Purity (JP method)	-
Related substances (1) (JP method)	to pass test
Related substances (2) (JP method)	to pass test
Assay (JP method)	min.98.0%
* Solubility in methanol	to pass test
* Loss on drying over silica gel	max.1.0%

\*:Additional test performed by Wako