

SPECIFICATION

Rhein

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	Yellow to reddish yellow powder
Identification (JP method)	to pass test
Absorbance E1%1cm (257nm) (JP method)	678~720
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
* Assay (qNMR) (JP method)	min.95.0%
** Solubility in acetone	to pass test
** Loss on drying over silica gel in vacuum	max.2.0%

*:The specification is originally set by Wako

**:Additional test performed by Wako