

Code No.19-0029

## SPECIFICATION

## Selenium

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Dark red ~ black, powder
Assay (Se)	min.99.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Nitrogen compounds (as N)	max.0.003%(mass/mass)
Sulfur compounds (as S)	max.0.02%(mass/mass)
Copper (Cu)	max.0.005%(mass/mass)
Lead (Pb)	max.0.005%(mass/mass)
Tellurium (Te)	max.0.03%(mass/mass)
Iron (Fe)	max.0.005%(mass/mass)

JIS K 8598:2018

\*:Additional test performed by Wako but not required by JIS

This product is powder