

SPECIFICATION

Silver Nitrate
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (AgNO ₃) Solubility in water Acidity Chloride (Cl) Sulfate (SO ₄) Substances not precipitated by hydrochloric acid (as sulfate) Copper (Cu) Lead (Pb) Iron (Fe)	White, crystals ~ crystalline powder min.99.8%(mass/mass) to pass test to pass test max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.01%(mass/mass) max.2ppm(mass/mass) max.2ppm(mass/mass) max.2ppm(mass/mass)

JIS K 8550:2021(*:Additional test performed by Wako but not required by JIS)