

SPECIFICATION

Sodium Acetate Trihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH ₃ COONa·3H ₂ O) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Phosphate (PO ₄) Sulfate (SO ₄) Potassium (K) Copper (Cu) Calcium (Ca) Aluminium (Al) Lead (Pb) Iron (Fe) Substances reducing permanganate (as O)	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 7.5~9.0 max.5ppm(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass) max.0.01%(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass)

JIS K 8371:2006*:Additional test performed by Wako but not required by JIS