

SPECIFICATION

Sodium Acetate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH ₃ COONa) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Phosphate (PO ₄) Sulfate (SO ₄) Potassium (K) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe) Substances reducing permanganate	White, crystals ~ crystalline powder min.98.5%(mass/mass) to pass test 7.5~9.0 max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.003%(mass/mass) max.0.02%(mass/mass) max.5ppm(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass)

JIS K 8372:2013

(*:Additional test performed by Wako but not required by JIS)