

Code No.19-0126

## SPECIFICATION

### Sodium L(+)-Aspartate Monohydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystalline powder
Solubility in water	to pass test
Specific rotation $[\alpha]_{D20}(c=2, \text{dil.HCl})$	+18.0~+22.0°
Loss on drying over silica gel in vacuum	max.0.3%
Thin layer chromatography test	to pass test
Assay (after drying)	min.97.0%