

SPECIFICATION

Disodium Glycerophosphate 5.5-Hydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals ~ powder
Solubility in water	to pass test
Chloride (Cl)	max.0.005%
Phosphate (PO ₄)	max.0.05%
Sulfate (SO ₄)	max.0.01%
Heavy metals (as Pb)	max.0.002%
Iron (Fe)	max.0.001%
Assay	97.0~103.0%