

Code No.19-0248

# SPECIFICATION

## Disodium Molybdate(VI) Dihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ powder
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	6.0~10.0
Chloride (Cl)	max.0.02%
Nitrate	to pass test (NO <sub>3</sub> :約0.005%以下)
Phosphate (PO <sub>4</sub> )	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.02%
Lead (Pb)	max.0.002%
Iron (Fe)	max.0.004%
Ammonium (NH <sub>4</sub> )	max.0.02%
Assay	min.98.5%