

SPECIFICATION

Sodium Nitrite
Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--|
| * Appearance Assay (NaNO ₂) Solubility in water Loss on drying at 105°C Chloride (Cl) Sulfate (SO ₄) Potassium (K) Magnesium (Mg) Calcium (Ca) Lead (Pb) Iron (Fe) | White ~ slightly yellow, crystals ~ crystalline powder min.98.5%(mass/mass) to pass test max.0.5%(mass/mass) max.0.003%(mass/mass) max.0.002%(mass/mass) max.0.003%(mass/mass) max.0.002%(mass/mass) max.0.002%(mass/mass) max.2ppm(mass/mass) max.2ppm(mass/mass) |

JIS K 8019:2010*:Additional test performed by Wako but not required by JIS