

Code No.19-0360

## SPECIFICATION

## Sodium Thiosulfate

Extra Pure

REQUIREMENT	SPECIFICATION
* Appearance	White crystalline powder
Assay (Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> )	min.95.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	6.0~8.5
Sulfate and sulfite (as SO <sub>4</sub> )	max.0.4%(mass/mass)
Sulfide (S)	max.0.003%(mass/mass)
Magnesium (Mg)	max.0.002%(mass/mass)
Calcium (Ca)	max.0.004%(mass/mass)

JIS K 8638:2011(\*:Additional test performed by Wako but not required by JIS)