

Code No.19-0447

SPECIFICATION

Amidosulfuric Acid

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (HOSO ₂ NH ₂)	min.99.5%(mass/mass)
Solubility in water	to pass test
Residue after ignition (as sulfate)	max.0.005%(mass/mass)
Chloride (Cl)	max.5ppm(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Lead (Pb)	max.5ppm(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8587 : 2018

*:Additional test performed by Wako but not required by JIS