

Code No.19-0717

## SPECIFICATION

pH3.25 Sodium Citrate Buffer Solution (10 times concentrated)  
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* pH at 25°C	3.21~3.29
Concentration of Cl (Cl <sup>-</sup> )	1.410~1.450mol/l
Concentration of Na (Na <sup>+</sup> )	1.85~2.05mol/l
2,2'-Thiodiethanol	4.8~5.2vol%
* Ninhydrin reactive substances	to pass test

\*:Specification of 10 times diluted solution