

Code No.19-0719

SPECIFICATION

pH5.28 Sodium Citrate Buffer Solution (10 times concentrated)
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* pH at 25°C	5.24~5.32
Concentration of Cl (Cl-)	0.942~0.982mol/l
Concentration of Na (Na+)	3.4~3.6mol/l
* Ninhydrin reactive substances	to pass test

*:Specification of 10 times diluted solution